

Speed-of-Sound Mapping for Pulse-Echo Ultrasound Raw Data using Linked-Autoencoders

Farnaz Khun Jush¹, Peter M. Dueppenbecker², Andreas Maier¹

¹ Pattern Recognition Lab, Friedrich-Alexander-University, Erlangen, Germany

² Technology Excellence, Siemens Healthcare GmbH, Erlangen, Germany

Workshop on Machine Learning for Multimodal Healthcare Data (ML4MHD 2023), 29th of July, Honolulu



Motivation

Breast cancer is the leading cause of cancer-related death for women worldwide


Gold-Standard: Xray mammography

Disadvantage: reduced sensitivity for dense breasts

Mortality rates among cancers¹



(a) Digital mammogram image fails to reveal the mass (triangle denotes abnormality)²




(b) Handheld ultrasound easily depicts a mass due to invasive ductal carcinoma

¹Sung, H., Ferlay, J., Siegel, R.L., Laversanne, M., Soerjomataram, I., Jemal, A. and Bray, F., 2021. Global cancer statistics 2020: GLOBOCAN estimates of incidence and mortality worldwide for 36 cancers in 185 countries. CA: a cancer journal for clinicians.

²Thigpen, D., Kappler, A. and Brem, R., 2018. The role of ultrasound in screening dense breasts—a review of the literature and practical solutions for implementation. Diagnostics, 8(1), p.20.

Ultrasound Shortcomings

- Dependent on operator's expertise
- Qualitative
- Similar tissue interpretations



(a) Breast
Carcinoma,
Malignant and
cancerous³

(b) Fibroadenoma,
Benign and non-
cancerous


³Ultrasound of the Breast, Robin Smithuis, Lidy Wijers and Indra Dennert Alrijne, hospital in Leiderdorp - the Netherlands , <https://radiologyassistant.nl/breast/ultrasound/ultrasound-of-the-breast>

Motivation

Beyond B-mode imaging


Quantitative Information

- Density
- Speed-of-sound
- Attenuation



Quantitative Information

- Density
- Speed-of-sound
- Attenuation



⁴Li, C., Duric, N., Littrup, P. and Huang, L., 2009. In vivo breast sound-speed imaging with ultrasound tomography. *Ultrasound in medicine & biology*, 35(10), pp.1615-1628.




Friedrich-Alexander-Universität
Technische Fakultät





Prior Works

- **Data Acquisition:** Transducer encodes an intermediate representation of the object under examination in the sensor domain.
- **Reconstruction:** Inversion of the corresponding encoding function:




End-to-end Networks (domain transfer) trained on simulated data :


Take Raw data (RF data as input) and returns speed-of-sound map in the output



(a) Raw data




(b) Network architecture consists of encoding and decoding paths



(c) Speed-of-sound ground truth

Prior works

Stability Challenges






Methods

Autoencoders:



- I. Encode input data to an intermediate representation
- II. From the intermediate representation reconstructs input in their outputs




Training

Optimizer: Adam

L2 loss: $X \rightarrow X'$, $Y \rightarrow Y'$, $M \rightarrow N$








Results

Results


Over one frame, same field of view




Results

Over one frame, same field of view


Mean SoS value inside the inclusion




Mean SoS value in the background




Stability comparison: AutoSpeed vs. En-De-Net



(a)



(b)



(c)


(a): 20 frames of same field of view from CIRS 073 breast mimicking phantom of size 3.8x3.8 cm in lateral and depth
(b): Corresponding reconstructed SoS maps using En-De-Net, (c): Corresponding reconstructed SoS maps using AutoSpeed
AutoSpeed shows more stability in comparison to the end-to-end trained network

Results

Multiple frames, dense inclusion in different locations

Networks Comparison:

	Background (m/s)	Inclusion (m/s)
Expected Value	1520 ± 10	1580 ± 20
AutoSpeed	1535 ± 6	1561 ± 11
En-De-Net	1527 ± 19	1545 ± 45



Conclusion

- We showed that SoS mapping is possible by employing a Linked autoencoder setup
- We showed that the linked autoencoder approach is more stable compared to the end-to-end mapping solution previously being used
- Even though the network is trained using only simulated data, on the measured data setup, the predicted SoS maps are close to the expected range

Outlook

- More research is required to transfer and prove the efficiency of such methods in clinical setups



Friedrich-Alexander-Universität
Technische Fakultät

SIEMENS
Healthineers



Thank you for your Attention!